Malate Assay Kit (Colorimetric)

LS-K176-100 (100 Tests) • Store at -20°C



Introduction

L-malate, or L-Malic acid, is a dicarboxylic acid that is made by all living organisms and plays an important role in the Calvin and Krebs Cycle. It is a source of CO₂ for the Calvin cycle in plants and is also an intermediate that forms from fumarate in the Krebs Cycle. Malate is frequently used in food and beverage industries as an additive in products such as wine, beer, candies, etc.

This L-malate assay kit is based on malate dehydrogenase catalyzed oxidation of malate in which the formed NADH reduces a formazan (MTT) reagent. The intensity of the product color, measured at 565 nm is proportional to the malate concentration in the sample.

Key Features

- Fast and sensitive. Use of 20 μL sample. Linear detection range 0.02 to 2 mM L-malate in 96-well plate assay.
- Convenient. The procedure involves adding a single working reagent, and reading the optical density at time 15 minutes. Room temperature assay. No 37°C heater is needed.
- High-throughput. Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.

Applications

Direct Assays: malate in food, juice, beverage and other agricultural products.

Components

	K176-100	
Component	100 Tests	
Assay Buffer	10 mL	
NAD/MTT	1 mL	
Enzyme A	120 μL	
Enzyme B	120 μL	
Standard (20mM L-Malate)	1.0 mL	

Materials Not Supplied

Pipetting devices and accessories (e.g. multi-channel pipettor), clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760), centrifuge tubes and plate.

Storage

The kit is shipped on ice. Store all components at -20°C upon receiving. Shelf life: 6 months after receipt.

Malate Assay Kit (Colorimetric)

LS-K176-100 (100 Tests) • Store at -20°C



Assay Procedure

Sample Preparation

Clear and slightly colored samples can be assayed directly.

Solid samples (food, fruits, etc.) can be homogenized in water followed by filtration or centrifugation (e.g. 5 min 14,000 rpm). Beverage samples can be assayed directly. Prior to assay, check the pH of the sample. If the pH is not between 7 and 8, adjust the sample pH to 7-8 with NaOH or HCl. Samples containing carbon dioxide should be degassed by gentle stirring prior assay. No dilution necessary in general.

It is prudent to test several dilutions to determine an optimal dilution factor *n*.

All samples can be stored at -20 to -80°C for at least one month.

Procedure

1. Standards and Samples. Equilibrate all components to room temperature. Prepare 500 μ L 2.0 mM L-Malate Premix by mixing 50 μ L 20 mM Standard and 450 μ L distilled water. Dilute standard as follows.

No	Premix + H ₂ O	Vol (µL)	L-Malate (mM)
1	100 μL + 0 μL	100	2.0
2	60 μL + 40 μL	100	1.2
3	30 μL + 70 μL	100	0.6
4	0 μL + 100 μL	100	0

Transfer 20 μ L standards into separate wells of a clear, flat-bottom 96-well plate. Transfer 20 μ L of each sample into two separate wells, one serving as a sample blank well (OD_{BLANK}) and one as a sample well (OD_{SAMPLE}).

2. Prepare sufficient Working Reagent (WR) by mixing for each standard and sample well, $74~\mu$ L Assay Buffer, $1~\mu$ L Enzyme A, $1~\mu$ L Enzyme B and $8~\mu$ L NAD/MTT. Prepare blank Working Reagent (BWR) by mixing for each sample blank well, $74~\mu$ L Assay Buffer, $1~\mu$ L Enzyme B and $8~\mu$ L NAD/MTT (i.e. no Enzyme A). Fresh reconstitution of the WRs is recommended.

Add 80 μ L WR to the four Standards and the Sample Wells. Add 80 μ L BWR to the Sample Blank Wells. Tap plate to mix briefly and thoroughly. Incubate 15 minutes at room temperature.

3. Read optical density at 565nm (520-600nm).

Calculations

Subtract the blank value (#4) from the standard values and plot the Δ OD against standard concentrations. Determine the slope and calculate the L-malate concentration of Sample,

[L-Malate] =
$$\frac{OD_{SAMPLE} - OD_{BLANK}}{Slope (\mu M^{-1})} \times n (\mu M)$$

 $\mathsf{OD}_{\mathsf{SAMPLE}}$ and $\mathsf{OD}_{\mathsf{BLANK}}$ are optical density readings of the Sample and Sample Blank, respectively. n is the sample dilution factor.

Note: if the sample OD value is higher than OD for 2 mM L-malate standard, dilute sample in water and repeat the assay. Multiply the results by the dilution factor.


Conversions: 1 mM L-malate equals 13.3 mg/dL, 0.018% or 133 ppm.

Malate Assay Kit (Colorimetric)

LS-K176-100 (100 Tests) • Store at -20°C



Sample Data



Standard Curve in 96-well plate assay in water.

Version: V.08.09.2018