

Xanthine Assay Kit (Colorimetric/Fluorometric)

LS-K200-100 (100 Tests) • Store at -20°C



Introduction

Xanthine is a purine base that can be found in most animal tissues and fluids. It is a product in the purine degradation pathway, produced by guanine deaminase from guanine, and by xanthine oxidoreductase from hypoxanthine. Xanthine is degraded to uric acid by xanthine oxidase. Clinically, xanthine and its derivatives act on sleep-inducing adenosine receptors as antagonists. Simple, direct and high-throughput assays for measuring xanthine find wide applications in research and drug discovery. This xanthine assay kit uses a single Working Reagent that combines the xanthine oxidase reaction and color reaction in one step. The change in color intensity of the reaction product at 570 nm or fluorescence intensity at $\lambda_{ex/em} = 530/585$ nm is directly proportional to xanthine concentration in the sample.

Key Features

- Sensitive and accurate. Use as little as 10 μ L samples. Linear detection range in 96-well plate for 30 minute incubation: 0.01 to 2.5 mM xanthine for colorimetric assays and 3 to 250 μ M for fluorometric assays.
- Simple and convenient. The procedure involves addition of a single working reagent and incubation for 30 min at room temperature.
- Fast and high-throughput. Assays using 96-well plates and liquid handling system could allow simultaneous processing tens of thousands of samples per day.

Applications

- Direct Assays: xanthine concentration in cell lysate, serum, and other biological samples.
- Drug Discovery/Pharmacology: effects of drugs on xanthine (purine) metabolism.

Components

| Component | K200-100 |
|--------------------------|-------------|
| | 100 Tests |
| Assay Buffer | 10 mL |
| HRP Enzyme | 120 μ L |
| Standard (2 mM Xanthine) | 1 mL |
| Dye Reagent | 120 μ L |
| XO Enzyme | 100 μ L |

Materials Not Supplied

Pipetting devices, centrifuge tubes, clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760), black 96-well plates (e.g. Greiner Bio-One, cat# 655900) and plate reader capable of either measuring absorbance between 550-585 nm or fluorescence intensity at $\lambda_{ex/em} = 530/585$ nm.

Storage

The kit is shipped on ice. Store all reagents at -20°C. Shelf life: 6 months after receipt.

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121
www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

Xanthine Assay Kit (Colorimetric/Fluorometric)

LS-K200-100 (100 Tests) • Store at -20°C



Assay Procedure

Colorimetric Procedure

Samples can be analyzed immediately after collection, or stored in aliquots at -20°C. Avoid repeated freeze-thaw cycles. If particulates are present, centrifuge sample and use clear supernatant for assay.

1. Equilibrate all components to room temperature. During experiment, keep thawed Enzyme in a refrigerator or on ice.
2. Standard Curve. Prepare standards as shown in the Table below.

| No | Standard + H ₂ O | Vol (μL) | Xanthine (mM) |
|----|-----------------------------|----------|---------------|
| 1 | 100 μL + 0 μL | 100 | 2 |
| 2 | 60 μL + 40 μL | 100 | 1.2 |
| 3 | 30 μL + 70 μL | 100 | 0.6 |
| 4 | 0 μL + 100 μL | 100 | 0 |

Transfer 10 μL standards and samples into separate wells.

3. Working Reagent. Prepare bulk working reagent by mixing 90 μL Assay Buffer, 1 μL XO Enzyme, 1 μL HRP Enzyme (vortex briefly before pipetting), and 1 μL Dye Reagent per reaction well in a clean tube. Transfer 90 μL Working Reagent into each reaction well. Tap plate to mix.
4. Incubate 30 min at room temperature, and then read optical density at 570 nm (550-585 nm) (OD₃₀).

Fluorometric Procedure

For fluorometric assays, the linear detection range is 3 to 250 μM xanthine. Dilute the standards from Colorimetric Procedure 10x with dH₂O to obtain standards at 200, 120, 60 and 0 μM Xanthine.

Transfer 10 μL standards and 10 μL samples into separate wells of a black 96-well plate.

Add 90 μL Working Reagent (see Colorimetric Procedure), tap plate to mix.

Incubate 30 min at room temperature, and then read fluorescence at λ_{ex/em} = 530/585 nm (F₃₀).

Calculations

Subtract blank OD₃₀ or F₃₀ (water, #4) from all standards and samples OD₃₀ or F₃₀ values and plot the ΔOD or ΔF against standard concentrations. Calculate the concentration using the equation below:

$$[\text{Xanthine}] = \frac{R_{\text{SAMPLE}} - R_{\text{BLANK}}}{\text{Slope } (\mu\text{M}^{-1})} \times n \quad (\mu\text{M})$$

Where R_{Sample} and R_{Blank} are the optical density or fluorescent values of the sample and blank, respectively. Slope is the slope of the standard curve and *n* is the dilution factor.

Notes: If the calculated sample xanthine concentration is higher than 2 mM in colorimetric assay or 200 μM in fluorometric assay, dilute sample in water and repeat the assay. Multiply result by the dilution factor (*n*).

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121



www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

Xanthine Assay Kit (Colorimetric/Fluorometric)

LS-K200-100 (100 Tests) • Store at -20°C



Sample Data



Version: V.08.09.2018

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com