

# NAD<sup>+</sup>/NADH Assay Kit (Fluorometric)

LS-K263-100 (100 Tests) • Store at -20°C



---

## Introduction

Pyridine nucleotides play an important role in metabolism and, thus, there is continual interest in monitoring their concentration levels. Quantitative determination of NAD<sup>+</sup>/NADH has applications in research pertaining to energy transformation and redox state of cells or tissue. This NAD<sup>+</sup>/NADH assay kit is based on a lactate dehydrogenase cycling reaction, in which the formed NADH reduces a probe into a highly fluorescent product. The fluorescence intensity of this product, measured at  $\lambda_{ex/em} = 530/585$  nm, is proportional to the NAD<sup>+</sup>/NADH concentration in the sample. This assay is highly specific for NAD<sup>+</sup>/NADH with minimal interference (<1%) by NADP<sup>+</sup>/NADPH and is a convenient Method to measure NAD, NADH and their ratio.

## Key Features

- Sensitive and accurate. Detection limit of 0.02  $\mu$ M and linearity up to 1  $\mu$ M NAD<sup>+</sup>/NADH in 96-well plate assay.
- Convenient. The procedure involves adding a single working reagent, and reading the fluorescence at time zero and 10 min.
- High-throughput. Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.

## Applications

- Direct Assays: NAD<sup>+</sup>/NADH concentrations and ratios in cell or tissue extracts.

## Components

Component	K263-100
	100 Tests
Assay Buffer	10 mL
Lactate	1.5 mL
Probe	750 $\mu$ L
NAD/NADH Extraction Buffers	12 mL each
Enzyme A	120 $\mu$ L
Enzyme B	120 $\mu$ L
NAD Standard	0.5 mL

## Materials Not Supplied

Pipetting (multi-channel) devices. Black, flat-bottom 96-well plates and fluorescent plate reader capable of reading at  $\lambda_{ex/em} = 530/585$  nm.

## Storage

The kit is shipped on ice. Store all reagents at -20°C. Shelf life: 6 months after receipt.

---

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# NAD<sup>+</sup>/NADH Assay Kit (Fluorometric)

LS-K263-100 (100 Tests) • Store at -20°C



## Assay Procedure

Note: This kit can also be used directly on cells cultured in 96 well plates.

### Procedure

1. Sample Preparation. For tissues weigh ~20 mg tissue for each sample, wash with cold PBS. For cell samples, wash cells with cold PBS and pellet  $\sim 10^5$  cells for each sample. Homogenize samples (either tissue or cells) in a 1.5 mL Eppendorf tube with either 100  $\mu$ L NAD extraction buffer for NAD determination or 100  $\mu$ L NADH extraction buffer for NADH determination. Heat extracts at 60°C for 5 min and then add 20  $\mu$ L Assay Buffer and 100  $\mu$ L of the opposite extraction buffer to neutralize the extracts. Briefly vortex and spin the samples down at 14,000 x g for 5 min. Use supernatant for NAD/NADH assays. Determination of both NAD and NADH concentrations requires extractions from two separate samples.
2. Calibration Curve. Prepare 5000  $\mu$ L 1  $\mu$ M NAD Premix by mixing 5  $\mu$ L 1 mM Standard and 4995  $\mu$ L distilled water. Dilute standard as follows.

No	Premix + H <sub>2</sub> O	NAD ( $\mu$ M)
1	100 $\mu$ L + 0 $\mu$ L	1.0
2	60 $\mu$ L + 40 $\mu$ L	0.6
3	30 $\mu$ L + 70 $\mu$ L	0.3
4	0 $\mu$ L + 100 $\mu$ L	0

Transfer 50  $\mu$ L standards into wells of a black flat-bottom 96-well plate.

3. Samples. Add 50  $\mu$ L of each sample in separate wells.
4. Reagent Preparation. For each reaction well, prepare Working Reagent by mixing 40  $\mu$ L Assay Buffer, 1  $\mu$ L Enzyme A, 1  $\mu$ L Enzyme B, 10  $\mu$ L Lactate and 5  $\mu$ L Probe. Fresh reconstitution is recommended.
5. Reaction. Add 50  $\mu$ L Working Reagent per well quickly. Tap plate to mix.
6. Read fluorescence at  $\lambda_{ex/em} = 530/585$  nm for time “zero” ( $F_0$ ) and  $F_{10}$  after a 10-min incubation at room temperature. Protect plate from light during this incubation.

### Calculations

First compute the  $\Delta F$  for each standard and sample by subtracting  $F_0$  from  $F_{10}$ . Plot the standard  $\Delta F$ 's and determine the slope. The NAD(H) concentration of the sample is computed as follows:

$$[\text{NAD(H)}] = \frac{\Delta F_{\text{SAMPLE}} - \Delta F_{\text{BLANK}}}{\text{Slope } (\mu\text{M}^{-1})} \times n \quad (\mu\text{M})$$

where  $\Delta F_{\text{SAMPLE}}$  and  $\Delta F_{\text{BLANK}}$  are the change in fluorescence intensity values of the Sample and Blank (STD 4) respectively. Slope is the slope of the standard curve and  $n$  is the dilution factor (if necessary).

Note: If the sample  $\Delta F$  values are higher than the  $\Delta F$  value for the 1  $\mu$ M standard, dilute sample in distilled water and repeat this assay. Multiply the results by the dilution factor.

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# NAD<sup>+</sup>/NADH Assay Kit (Fluorometric)


LS-K263-100 (100 Tests) • Store at -20°C

---

## General Considerations

1. At these concentrations, the standard curves for NAD and NADH are identical. Since NADH in solution is unstable, we provide only NAD as the standard.
2. This assay is based on an enzyme-catalyzed kinetic reaction. Addition of Working Reagent should be quick and mixing should be brief but thorough. Use of multi-channel pipettor is recommended.
3. The following substances interfere and should be avoided in sample preparation. EDTA (>0.5 mM), ascorbic acid, SDS (>0.2%), sodium azide, NP-40 (>1%) and Tween-20 (>1%).
4. For samples containing higher than 100  $\mu$ M pyruvate, we recommend using an internal standard.

## Sample Data



Version: V.08.09.2018

---

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)