Monoamine Oxidase (MAO) Assay Kit (Fluorometric)

LISBio LifeSpan BioSciences, Inc.

LS-K332-100 (100 Tests) • Store at -20°C

Introduction

Monoamine oxidases (MAO, EC 1.4.3.4) are a family of mitochondrial enzymes that catalyze the oxidative deamination of monoamines. MAO dysfunction is thought to be responsible for a number of neurological disorders. Unusually high or low levels of MAOs in the body have been associated with depression, schizophrenia, substance abuse, attention deficit disorder, migraines, and irregular sexual maturation. MAO inhibitors are one of the major classes of drug prescribed for the treatment of depression. This MAO Assay Kit provides a convenient fluorometric means to measure MAO enzyme activity. In the assay, MAO reacts with p-tyramine, a substrate for both MAO-A and MAO-B, resulting in the formation of H_2O_2 , which is determined by a fluorometric Method ($\lambda_{ex/em} = 530/585$ nm). The assay is simple, sensitive, stable and high-throughput adaptable.

Key Features

- Safe. Non-radioactive assay.
- Sensitive and accurate. As low as 0.01 U/L MAO activity can be quantified.
- Homogeneous and convenient. "Mix-incubate-measure" type assay. No wash and reagent transfer steps are involved.
- Robust and amenable to HTS: can be readily automated on HTS liquid handling systems for processing thousands of samples per day.

Applications

- MAO-A/B activity determination in biological samples.
- Evaluation and screening for MAO inhibitors.

Components

	K332-100
Component	100 Tests
Assay Buffer (pH 7.4)	12 mL
Pargyline (20 mM)	50 μL
Clorgyline (20 mM)	50 μL
Hydrogen Peroxide (3% H ₂ O ₂)	100 μL
<i>p</i> -Tyramine	120 μL
HRP Enzyme	120 μL
Dye Reagent	120 μL

Materials Not Supplied

Pipetting devices, centrifuge tubes, black flat-bottom 96-well plate (e.g. Corning Costar).

Storage

The kit is shipped on ice. Store all components at -20°C. This product is shipped on ice. Shelf life: 6 months after receipt.

Monoamine Oxidase (MAO) Assay Kit (Fluorometric)

LISBio
LifeSpan BioSciences, Inc.

LS-K332-100 (100 Tests) • Store at -20°C

Assay Procedure

Notes

- (1). Thiols (β -mercaptoethanol, dithioerythritol etc) at > 10 μ M interfere with this assay and should be avoided in sample preparation.
- (2). Samples should be free of particle or precipitates. MAO can be extracted from a tissue by homogenization and differential centrifugation, e.g. Biochem. J. (1968) 108: 95. Store sample at -80°C.
- (3). Prior to assay, concentrations of protein, inhibitor, substrate and incubation time may need to be established for a given sample.

Procedure

Use black flat-bottom plates. Prior to assay, bring all components to room temperature, briefly centrifuge tubes before opening.

Dilute the 20 mM inhibitors with H_2O to 10 μ M (e.g. mix 5 μ L 20 mM inhibitor with 10 mL H_2O).

- 1. To determine MAO-A activity, use 1 mM p-tyramine substrate and include a control with 0.5 μ M MAO-A inhibitor clorgyline.
 - Samples: dilute sample in Assay Buffer. Transfer 45 μ L of each sample into two separate wells. Add 5 μ L H₂O (SAMPLE) and 5 μ L 10 μ M clorgyline (CONTROL). Mix and incubate for 10 min at room temperature for the inhibitor to block MAO-A activity.
- 2. Calibrator. Mix 5 μ L H₂O₂ with 1400 μ L H₂O. Further dilute 5 μ L of the resulting H₂O₂ in 780 μ L H₂O to give 20 μ M H₂O₂. Dilute calibrator with H₂O to give 20, 10, 5 and 0 μ M H₂O₂.
 - Transfer 50 µL calibrators into separate wells of the assay plate.
- 3. Prepare enough Working Reagent for all sample and calibrator wells. For each well, mix: $50~\mu$ L Assay Buffer, $1~\mu$ L ptyramine, $1~\mu$ L Dye Reagent and $1~\mu$ L HRP Enzyme. Transfer $50~\mu$ L Working Reagent to all wells. Briefly tap plate to mix.
- 4. Incubate for 20 min in the dark. Read fluorescence intensity at λ_{exc} = 530 nm and λ_{em} = 585 nm.

To measure MAO-B activity, use 1 mM p-tyramine and include a control with 0.5 μ M pargyline (MAO-B inhibitor). Procedure is the same as for MAO-A determination.

To screen for MAO inhibitors or characterize inhibitor potency (IC₅₀), mix 5 μ L inhibitor with 45 μ L sample and incubate for at least 10 min to allow the inhibitor to interact with the enzyme, prior to adding the Working Reagent.

Calculations

Plot H₂O₂ calibration curve and determine its Slope (μM⁻¹). MAO enzyme activity in the sample is calculated as

MAO Activity =
$$\frac{RFU_{SAMPLE} - RFU_{CONTROL}}{Slope \times t}$$
 (U/L)

where RFU_{SAMPLE} and RFU_{CONTROL} are the measured fluorescence values of the sample and sample control (i.e., in the presence of the respective inhibitor pargyline or clorgyline). *t* is the incubation time (20 min).


Unit definition: one unit of MAO catalyzes the formation of 1 μmole H₂O₂ per min under the assay conditions.

Monoamine Oxidase (MAO) Assay Kit (Fluorometric)



LS-K332-100 (100 Tests) • Store at -20°C

Sample Data



Version: V.08.09.2018